

Suction Machine

Portable Type

— 9E-A



Oil-free piston pump



Stable working performance



Anti over-flow technology



Large pump rate



Silence and low noise



Portable and home use

Products Features

- Heavy duty oil-free piston pump
- Anti over-flow technology and large pump rate
- Silence and stable working performance
- 1000 ml polycarbonate bottle shatter-proof and washable
- Water proof switch and metal pump nozzle
- Portable design suitable for home and clinic use

Specifications

Input Voltage	AC220V± 10%
Input Frequency	50Hz± 2%
Average Power Consumption	90VA
Max Negative Pressure Value	≥ 0.075MPa
Pumping Rate	18L/min
Continuous Working Time	≥ 30 min
Rest Time	≤ 30 min
Weight	3.9 kg
Sound Level	≤ 53dB
Dimensions	35*26*29cm
Jug Volume	1000 ml
Package	2pcs/carton 54.5*38*32 cm

Accessories



Liquid storage bottle



Bacterial filter



Medical silicon tube



Dual piston pump

