

HEMOCARE OXYGEN CONCENTRATOR

MODEL: JAY-3AW/JAY-5AW



Low noise,
takes up less space



USA imported humidifier,
6PSI pressure release



Innovative big
LCD display

Optional

1. Nebulizer(Atomization)
2. Low Purity alarm
3. Low purity alarm with low flow alarm
4. High temperature alarm
5. Pulse oximeter
6. Maintenance alarm
7. Reset Switch
8. Remote control monitoring system (WIFI/Bluetooth)

Specification

Flow Rate	0-3LPM / 0-5LPM
Purity	93% ($\pm 3\%$)
Outlet Pressure	0.04-0.07MPA(6-10PSI)
Sound Level	≤ 40 db
Power Consumption	≤ 300 w / ≤ 350 w
LCD Display	Switch times, operating pressure, present working time, accumulated working time, presetting time from 10mins to 40hours
Alarm	Power Failure alarm, High&low pressure alarm
Net Weight	16Kgs / 16 Kgs
Size	348×280×510 mm